

| Name of Material | : | Bottle – Glass |
|-------------------|---|----------------|
| Manufacturer | : | NA |
| Supplier | : | NA |
| Country of Origin | : | NA |
| Batch/Lot No. | : | NA |
| MFG Date | : | NA |
| EXP Date | : | NA |
| Challan No. | : | NA |
| QA Ref. No. | : | NA |
| QRN Ref. No. | : | NA |
| Qty. Supplied | : | NA |
| Date of Received | : | NA |
| Date of Testing | : | NA |
| Date of Reporting | : | NA |

| SN | Description | Specification | Results |
|----|----------------------------------|-----------------------|---------|
| 01 | Bottle Weight | As per design | |
| 02 | Bottle Volume | As per design | |
| 03 | Mouth Dia | As per design | |
| 04 | Neck Dia | As per design | |
| 05 | Body Dia | As per design | |
| 06 | Bottom Dia | As per design | |
| 07 | Bottle Height | As per design | |
| 08 | Bottle Thickness – Neck Side | As per design | |
| 09 | Bottle Thickness – Body Side | As per design | |
| 10 | Bottle Thickness – Bottom Side | As per design | |
| 11 | Bottle Color | As per design | |
| 12 | Bottle Material Distribution | Smooth Distribution | |
| 13 | Bottle Leakage | No Leakage | |
| 14 | Bottle Straight Standing Ability | Must Stand Straightly | |
| 16 | Air Gap in Bottle Body | Must be Absent | |
| 17 | Milky Wake in Bottle Body | Must be Absent | |
| 18 | Bottle Body Leakage | Must be Absent | |



| | |
|----------------|--|
| Remarks | |
|----------------|--|

| | | | | |
|-----------------|--|----------|--|----------|
| Decision | | Accepted | | Rejected |
|-----------------|--|----------|--|----------|

.....
Tested By

.....
Checked By

.....
Approved By



Other Related Information

| Name of Material | Bottle – Glass |
|--------------------------|--|
| Others Name | NA |
| Chemical Composition | NA |
| Molar Mass | NA |
| Density | NA |
| Melting Point | NA |
| Boiling Point | NA |
| Others Information | NA |
| Function in Food Process | <ul style="list-style-type: none">- Used in Juice Filling line- Used in CSD line- Used in oil bottling line- Used for other type of liquid filling line- Used in pharmaceuticals |